


Fisherbrand
QUALITY. RELIABILITY. VALUE.

Polypropylene & Polystyrene Plates

Labware and Life Science Products



Fisherbrand
QUALITY. RELIABILITY. VALUE.

Polypropylene & Polystyrene Plates


Labware and Life Science Products

Microplate readers, Multiskan™ GO

- Very fast plate measurements and a full spectrum of a sample in less than 10sec
- Incubation up to 45°C and variable shaking speeds to optimise sensitive applications such as enzyme kinetics and cellular assays
- A unique power save function for reduced energy consumption

Applications:

- Immunoassays
- Protein assays, enzymes assays
- Growth curve and hormone assays
- Endotoxins, cytotoxicity
- Food diagnostics
- HIV assays
- Hybridisation assays
- Mini sequencing assays
- Cell proliferation, cell adhesion
- Signal transduction



11131_2015

© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: (0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 2304 932-5 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500


 **Fisher
Scientific**
A Thermo Fisher Scientific Brand

 **Fisher
Scientific**

Polystyrene Plates

Fisherbrand Polystyrene Microplates are optimized for manual or robotic assays

- 96-well microplates are ideal for antibiotic screens, cell-based assays and screening compounds
- 384-well microplates conserve samples and reagents, providing cost savings and greater screening productivity
- 1536-well microplates are engineered for high throughput screening, allowing end users to screen four times as many samples in one plates



Fisherbrand Cat. No.	Description	Case Qty
11987944	Fisherbrand 96 Well Plate Round Bottom , Clear, Polystyrene, Non-Treated, Non-Sterile	80
11997944	Fisherbrand 96 Well Flat Bottom Plate, Clear, Polystyrene, Non-Treated, Non-Sterile	80
11947954	Fisherbrand 384 Well Flat Bottom Plate, Clear Polystyrene, Non-Treated, Non-Sterile	80
13525440	Fisherbrand 96 Well Flat Bottom Plate,White, Polystyrene, Non-Treated, Non-Sterile, w/o Lid	180
13535440	Fisherbrand 96 Well Flat Bottom Plate, Black, Polystyrene, Non-Treated, Non-Sterile, w/o Lid	180
13575480	Fisherbrand 96 Well Round Bottom Plate, Clear, Polystyrene, Non-Treated, Non-Sterile, w/o Lid	180
13545440	Fisherbrand 384 Well Flat Bottom Plate, White, Polystyrene, Non-Treated ,Non-Sterile, w/o lid	100
13555440	Fisherbrand 384 Well Flat Bottom Plate, Black, Polystyrene, Non-Treated, Non-Sterile, w/o Lid	100
13565440	Fisherbrand 384 Well Flat Bottom Plate, Clear Polystyrene, Non-Treated, Non-Sterile, w/o Lid	100
13575440	Fisherbrand 1536 Well Flat Bottom Plate, Clear Polystyrene, Non-Treated ,Non-Sterile, w/o Lid	100
13585440	Fisherbrand 1536 Well Flat Bottom Plate, Black Polystyrene, Non-Treated, Non-Sterile, w/o Lid	100
13595440	Fisherbrand 1536 Well Flat Bottom Plate, White Polystyrene, Non-Treated, Non-Sterile, w/o Lid	100

Accessories

Microplate lids protect samples from contamination and evaporation during assay processing, incubation or storage

Fisherbrand Cat. No.	Description	Case Qty
11927964	Fisherbrand Lid for 96/384 Well Plates, Clear, Polystyrene, Non-Sterile	180




Fisherbrand Polystyrene (PS) plates are optimized for robotic instruments, plate readers and liquid handling systems. The optical clarity of clear PS makes it the perfect material for many imaging applications. Background, autofluorescence and cross-talk are minimized with black or white PS plates.

Polypropylene Plates

Fisherbrand Polypropylene Plates are ideal for Storage Applications

- 96-well polypropylene (PP) microplates feature a shared-wall design which yields increased well volume
- 96-well deepwell microplates accommodate even larger working volumes - up to 1.0 mL or 2.0 mL storage volume per well
- 384-well microplates are ideal for small volume storage applications; conical wells provide optimal sample recovery
- 384-well deepwell microplates store more material and are ideal for compound library storage



Fisherbrand Polypropylene (PP) microplates are ideal for storage because of the low-binding characteristics of polypropylene. Proteins and DNA will not adhere to PP, allowing for complete sample recovery.

Fisherbrand Cat. No.	Description	Case Qty
11907954	Fisherbrand 96 Well Plate Round Bottom , Natural, Polypropylene, Non-Sterile	80
13505450	Fisherbrand 96 Well Round Bottom Plate, Natural, Polypropylene, Sterile	120
11917954	Fisherbrand 96 Well Plate Round Bottom , Red, Polypropylene, Non-Sterile	80
11927954	Fisherbrand 96 Well Plate Round Bottom, Yellow, Polypropylene, Non-Sterile	80
11937954	Fisherbrand 96 Well Plate Round Bottom , Blue, Polypropylene, Non-Sterile	80
13545450	Fisherbrand 96 DeepWell 1mL Plate Round Bottom, Natural, Polypropylene, Non-Sterile	50
13535450	Fisherbrand 96 DeepWell 1mL Plate Round Bottom, Natural, Polypropylene, Sterile	50
13555450	Fisherbrand 384 DeepWell 252ul, Natural, Polypropylene, Non-Sterile	60
13565450	Fisherbrand 384 Shallow Well Conical Bottom Plate, Natural, Polypropylene, Non-Sterile	100
13575450	Fisherbrand 384 Shallow Well Conical Bottom Plate, Natural, Polypropylene, Sterile	100
13595450	Fisherbrand 384 Shallow Well Conical Bottom Plate, Black, Polypropylene, Non-Sterile	100
13585450	Fisherbrand 384 Shallow Well Conical Bottom Plate, White, Polypropylene, Non-Sterile	100
11957954	Fisherbrand 384 Well Plate Round Bottom, Natural, Polypropylene, Non Sterile	120
11967954	Fisherbrand 384 Well Plate Round Bottom, Red, Polypropylene Non-Sterile	120
11977954	Fisherbrand 384 Well Plate Round Bottom, Yellow, Polypropylene Non-Sterile	120
11987954	Fisherbrand 384 Well Plate Round Bottom, Yellow, Polypropylene Non-Sterile	120