

Filter units, bottle top, nylon membrane, with receiver, Thermo Scientific Nalgene® MF75™ series



For sterilising media and other reagents used in tissue culture and extremely sensitive cell culture work.

The nylon membrane has extremely low levels of extractables (less than 0.5% by weight), is naturally hydrophilic, contains no wetting agents and is alcohol resistant.

The filter unit is colour coded; red collar indicates nylon membrane. Sterile, non-pyrogenic and non-cytotoxic.

Filter units - filter funnel, filter receiver and leak proof closure

Catalogue No	Legacy No	Alt. No	Upper capacity, mL	Receiver capacity, mL	Pore size, µm	Membrane diameter, mm	Pack qty
10762732	-	150-0020	150	150	0.2	50	12
10760802	-	153-0020	250	250	0.2	50	12
10790412	-	153-0045	250	250	0.45	50	12
10578471	-	151-4020	500	500	0.2	75	12
10211321	-	151-4045	500	500	0.45	75	12
10444261	-	154-0020	500	1,000	0.2	75	12
10169700	-	154-0045	500	1,000	0.45	75	12



Filter units, bottle top, SFCA membrane, Thermo Scientific Nalgene® MF75™ series



- Disposable tissue culture filters with 50mm, 75mm or 90mm surfactant free cellulose acetate (SFCA) membranes
- Securely screw onto 33mm or 45mm neck glass media bottles
- Sterile

Bottle top filter only - for 33mm neck bottles

Catalogue No	Legacy No	Alt. No	Capacity, mL	Pore size, µm	Membrane diameter, mm	Pack qty
10352911	-	290-3320	150	0.2	50	12
10027780	-	290-3345	150	0.45	50	12
10037790	-	291-3320	500	0.2	75	12
10529051	-	291-3345	500	0.45	75	12
10139800	-	292-3320	1,000	0.2	90	12

Bottle top filter only - for 45mm neck bottles

Catalogue No	Legacy No	Alt. No	Capacity, mL	Pore size, µm	Membrane diameter, mm	Pack qty
10139560	-	290-4520	150	0.2	50	12
10624281	-	290-4545	150	0.45	50	12
10201371	-	291-4520	500	0.2	75	12
10549241	-	291-4545	500	0.45	75	12
10790592	-	292-4520	1,000	0.2	90	12

